

Pacif-i™

Bluetooth Thermometer Pacifier

Instruction Booklet

BLUE MAESTRO'S STANDARD LICENSE AND DISCLAIMER

Copyright ©2017 Blue Maestro Limited. All Rights Reserved. The copyright in the software to which this document relates and in this document belong exclusively to Blue Maestro Limited and no permission to use, copy, modify, and distribute this software or this documentation is permitted unless with the express prior written authorisation of Blue Maestro Limited. The license to use the software to which this document relates is granted so long as such software resides in hardware provided by Blue Maestro Limited.

No claims, representations or warranties, whether expressed or implied, including but not limited to warranties of merchantability, fitness for a particular purpose, of title, or of non-infringement of third party rights, are made by Blue Maestro as to the safety, reliability, durability or performance of Blue Maestro's products. Blue Maestro is not responsible for any liabilities resulting from negligence, misuse, modification, or alterations to the product by the user. Furthermore, Blue Maestro accepts no liability whatsoever for the safety, reliability, durability or performance of any of its products. IN NO EVENT, REGARDLESS OF CAUSE, SHALL BLUE MAESTRO BE LIABLE FOR ANY INDIRECT, SPECIAL, INCIDENTAL, PUNITIVE OR CONSEQUENTIAL DAMAGES OF ANY KIND, WHETHER ARISING UNDER BREACH OF CONTRACT, TORT (INCLUDING NEGLIGENCE), STRICT LIABILITY OR OTHERWISE, AND WHETHER BASED ON THIS AGREEMENT OR OTHERWISE, EVEN IF ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

User assumes responsibility for correct operation of the product and any software associated with it. User assumes responsibility for determining the suitability of the product to the user's needs, for configuring and using the product to meet those needs, and for the proper placement/location of the product in the environment it is being used. User assumes responsibility for verifying and interpreting results obtained from product use.

BLUE MAESTRO RESERVES THE RIGHT TO REVISE ITS STANDARD LICENSE AND DISCLAIMER WITHOUT NOTICE AT ANY TIME.

Introduction

Congratulations on your purchase of Pacif-i™, the world's first Bluetooth Smart temperature thermometer for smartphones that monitors your baby's temperature while it is in their mouth.

When used properly, this high quality product will give you over 12 months of accurate, dependable service.

WARNING: Pacif-i is for oral use only and is intended for use with children aged 5 and under.

Compatibility

Pacif-i requires an iOS or Android smartphone or tablet that is Bluetooth 4.0 or above compatible and the runs iOS 7 or above or Android 5 or above operating systems. Please check the specifications of your smartphone or tablet before you begin.

While Pacif-i has been made to high standards of care and cleanliness, it is recommended that before you put it in your baby's mouth for the first time that you wash the teat of Pacif-i with warm soapy water, rinse and dry. Please see section below on cleaning and caring for your Pacif-i for more information.

Download the Apps

There are apps for both Android and iOS available on Google Play and the App Store respectively. Go to the relevant app store, in search type either "Blue Maestro" or "Pacif-i" to locate the correct app. Download it and install it on your Android or iOS device.

Ensure Bluetooth is turned on

Ensure you have Bluetooth switched on in your smartphone's or tablet's settings. Please note that Pacif-i does not bond with your smartphone or tablet, so it will never show in the list of Bluetooth devices in your settings.

Inspect prior to each use

Each time you give Pacif-i to your baby please check for any damage, such as broken or cracked plastic features or split or torn silicone nipple. Where any damage is noticed, please do not use and discard Pacif-i.

While Pacif-i has been made to high standards of care and cleanliness, it is recommended that before you put it in your baby's mouth for the first time that you wash the teat of Pacif-i in accordance with the instructions contained herein.

Cleaning your Pacif-i


Use warm soapy water or rubbing alcohol (70% isopropyl alcohol) to clean your pacifier. Do not wash the unit with any thinner or chemical solvent.

WARNING: Do not boil or place in a dishwasher and under no circumstances sterilise in a microwave. Pacif-i is rated IP65. It should not be left submerged in water nor should it be squeezed when submerged. This will likely damage your Pacif-i and will void the warrant.

Do not leave a soother in direct sunlight or near a source of heat, or leave in disinfectant (“sterilising solution”) for longer than recommended, as this may weaken the teat.

Turning your Pacif-i on and off

The front coloured surface of Pacif-i is a switch. Press to switch on and press and hold for three seconds to turn off. Turning your Pacif-i off is helpful for extending its battery life and is recommended at night and when it is not being used.



Accuracy and calibration

The accuracy of Pacif-i can be affected by the way your baby suckles on the pacifier, the prevailing ambient temperature and other factors such as whether your baby has recently eaten or drunk something hot or cold.

Also taking temperature in the mouth can provide a temperature up to 1°F/0.5°C less than other measurements, such as in-ear.

Although there is no singular ‘normal’ temperature reading, a temperature taken orally ranging between 97°F/36°C and 100°F/37.6°C is considered normal. Your child’s physician should be consulted to determine what temperature range the physician considers to be normal. The best method to determine your child’s normal temperature is to use the thermometer when the child is feeling well.

We suggest you calibrate Pacif-i by comparing it to another trusted clinical thermometer such as an infra-red in-ear thermometer and then adjusting for any difference using the calibration feature within the app. NOTE: An incorrectly calibrated Pacif-i can provide results not indicative of your baby’s actual temperature.

Pacif-i is a parental aid

Pacif-i is a parental aid and should not be used for medical purposes. It should not be used for clinical evaluations or for diagnosing illness. If your baby appears ill consider seeking medical advice even where your baby’s temperature appears normal when using Pacif-i. Also if Pacif-i indicates your baby has a temperature, use a clinical thermometer for diagnostic purposes before considering any treatment you may administer your baby.

Cleaning your Pacif-i

Use warm soapy water or rubbing alcohol (70% isopropyl alcohol) to clean your pacifier. Do not wash the unit with any thinner or chemical solvent. **WARNING:** Do not boil or place in a dishwasher and under no circumstances sterilise in a microwave. Pacif-i is rated IP65. It should not be left submerged in water nor should it be squeezed when submerged. This will likely damage your Pacif-i and will void the warranty.

For European Customers

Pacif-i complies with EN1400:2013 “Child use and care articles - Soothers for babies and young children - Safety requirements and test methods” and the relevant provisions of the Radio and Telecom-munication Terminal Equipment (R&TTE) Directive.

Pacif-i is not a medical device for the purposes of the Medical Devices Directive 93/42/EEC and should not be used for diagnostic purposes.



FOR YOUR CHILD’S SAFETY - WARNING!

Inspect carefully before each use. Pull the soother in all directions. Throw away at the first signs of damage or weakness. Only use dedicated soother holders tested to EN 12586. Never attach other ribbons or cords to a soother, your child may be strangled by them.

For US Customers

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference and
- (2) this device must accept any interference received, including interference that may cause undesired operation.

FCC ID: 2ACYPAC1

Pacif-i is exempt from FDA premarket notification requirements.

Support and Further Information

Please visit www.bluemaestro.com for support and further information on your Pacif-i.

Safe Disposal of Pacif-i

This product shall not be treated as household waste since it represents electronic equipment. Instead it shall be handed over to the applicable collection point for the recycling of electrical and electronic equipment. By ensuring this product is disposed of correctly, you will help prevent potential negative consequences for the environment and human health, which could otherwise be caused by inappropriate waste handling.